

Intuitive Operation. Seamless Workflow.




# Simplify Experiment Design—Leave More Time for Research

Olympus cellSens software simplifies your workflow:

- Intuitive user interface.
- User-specified configurations to fit your application.
- Seamless functionality from image capture through report creation.

Spend less time with your imaging software and make faster progress toward your research goals.

## Imaging




### Image Capture

Capture multicolor, time-lapse, and z-stack images with ease. Select the appropriate capture button, add relevant parameters, and click "Start." The Process Manager and Experiment Manager make multidimensional imaging easy.




## Processing



### Viewing and Processing

View your data in the layout that looks good to you. Take advantage of advanced image processing functions, such as stitching, extended focus, deconvolution, and unmixing to prepare your data for analysis.


## Analyzing



### Measurement and Analysis

cellSens features region of interest, phase analysis, and cell counting capabilities. Export raw measurement data to MS Excel or a cellSens workbook with a single click.

## Reporting



### Collaboration and Communication

Customizable database and reporting solutions enable active collaboration with colleagues and coworkers. These functions make it simple to manage, share, and distribute images and analysis.



### Microscopy Research with a Personal Touch

Olympus microscopes enable new imaging techniques and push the boundaries of resolution at all magnifications. Olympus cellSens software improves productivity with efficient acquisition workflow, image processing capabilities, and analytical strength. Centered around the needs of demanding customers, cellSens software is flexible, customizable, and designed to adapt as application requirements evolve.



# Reduce Clutter by Displaying Only the Tools and Windows You Need

## It's Time to Get Personal

Olympus has been at the forefront of microscopy for over 90 years, developing microscopes and systems supporting a broad range of applications. As a result, we understand that your research has individual requirements that require individual solutions. All cellSens software features an easy-to-use interface that you control and customize according to your needs.

## Dynamic Interface

Workflow efficiency requires careful definition of tasks and tools at every stage. cellSens software has a dynamic interface that helps provide the tools you need, which are clearly available at each step. Olympus has created a number of interface layouts, each developed with capabilities appropriate to the user's needs.

- **Acquisition Layout**

—Select between different acquisition processes and adjust the camera settings.

- **Processing Layout**

—Post-acquisition functions such as image processing, measuring, collecting data, and presenting the resulting statistics.

- **Count & Measure Layout**


—Manual and automated measurement and object counting.

- **Reporting Layout**

—Generate reports to document and share results.

- **Create Layout**

—Define customized layouts to suit any workflow.



## Camera Control Panel

Camera control and configuration is a central aspect of the cellSens workflow. Parameters such as exposure time and pixel binning are just two examples of typical camera settings available to optimize image quality. cellSens Entry and Standard include control and configuration for all Olympus digital cameras and microscopes. cellSens Dimension software includes precision control over research grade cameras when the most sensitive scientific grade detector is required.




## Dark Application Skin

The Dark Application Skin reduces ambient light from the computer monitor, enabling cellSens software users to better maintain darkened environments during imaging, with an interface design that remains highly visible and usable.



Dark skin



**Need Help? Online Assistance or Support is Just a Click Away**

**Display Only the Functions You Need on the Toolbar**

**Common Functions Grouped in a Single Tab**  
 All necessary functions are placed where you want. Layout tabs make it easy to select functions according to your workflow. For instance, display camera control features in your Acquisition layout and then remove them from view when you switch to the Processing layout.

**Display or Hide Windows as You Require, or Use Auto-hide for Clean Operation**


**Graphical Experiment Manager (GEM)**  
 The GEM enables users to design complex experiments by simply dragging and dropping icons onto the canvas.

**Create Flexible Workflow Toolbars for Repetitive Operations**  
 Design customized tool windows and create your own toolbars for easy access to the most useful functions. Add icons and text to optimize usability to enhance workflow efficiency.


**Functional Panels are Grouped in Tabs for Easier Access**

### Create your own Layout


Organize the tools and windows for the job at hand to create a functional layout that works best for you.



Full screen



Floating panels



Docked panels

# Solutions to Empower Your Research

## Scientific Research Need

## Our Answers

Quickly define complex experiments without programming

### Graphical Experiment Manager (GEM)


Design experiments by connecting various command icons. Image acquisition is available for up to 6 dimensions (XYZT multipoint).



Study cell dynamics with fluorescent markers

### Intensity Analysis


Visualize changes in intensity over time and save this information for later analysis. Ratio Analysis functions enable calibration, display, and analysis of live/stored data reflecting changes in the intensity ratio between two acquisition channels.



Measurement of cell count and confluency without stain

### Cell Confluency Measurement


Measure cell count and confluency using phase contrast images. Create cell growth curves with improved accuracy based on larger sample populations by automatically including multiple images within your analysis.



Improved image detail

### Deconvolution

Apply proprietary and highly efficient post-processing tools for both CCD and confocal imaging to enhance the ability to differentiate between objects. Choose between 2D and optional 3D blind deconvolution.




Kei Ito, Ph. D.  
Institute of Molecular and Cellular Biosciences, University of Tokyo

Automate repeating tasks

### Macro Manager

Use the macro manager to automate typical acquisition and data analysis workflows. Macro commands can be applied to multiple images simultaneously and can reduce the time required to complete complex imaging and image analysis.



## Clinical Research Needs

Simplification of image acquisition



## Our Answers

### Simple Layout

▶ MOVIE

The "Simple Layout" improves efficiency and workflows for all users from novice to expert. All image acquisition functions are easily accessible for intuitive operation. This enables even untrained users to obtain excellent results.





Image conferencing and consultation



### Conference mode

▶ MOVIE

Fit live and static imaging to the display for the best presentation. Various annotation and graphical annotation tools are at your fingertips.




Observe an entire large sample at once



### Panoramic Imaging

▶ MOVIE

Create seamless panoramic and accurately stitched images using a motorized stage. Real-time stitching mode produces wide area scans using a manual stage.




Speed up manual counting procedures



### Object Counting

Perform manual counts with user defined classes. Generate counts and proportions for each class at the click of a mouse.




Nuclei counting with variable thresholding



### Particle Analysis

Set threshold levels for nuclei counts, or calculate parameters such as tissue slice total area and area ratios.






# A Range of Easy-to-Use Functions Turn Your Findings into Compelling Presentations

## Image Capture

### Graphical Experiment Manager (GEM)

Dimension

GEM is a flexible drag-and-drop interface to build simple or complex experiments within cellSens software. Combine actions within specialized frames to dictate the order and priority of automation. Easily acquire multichannel, Z-stacks, and time-lapse imaging across one or more sample positions. Perform two-channel simultaneous imaging within GEM using the cellSens High-End Device solution. GEM permits users to interact with the system during long timelapse imaging without terminating the experiment.



Loop command

### Capture Multidimensional Images


Dimension

The Process Manager makes it easy to capture multicolored and multidimensional images with just a couple of clicks when imaging with a motorized stage.

Dimension + Multiposition

Standard + Multiposition

The optional multiposition solution is used to automatically capture multipoint and large area images.



Process Manager Setting

Multicolor images


### Real-time panoramic imaging

▶ MOVIE

Dimension

Standard + Manual Process

Create panoramic images in real-time by freely moving the manual stage with the Manual Process solution. Manual Process Control is available as an option for cellSens Standard software and included within cellSens Dimension software.



Stage control

Real-time stitching




Real-time stitched image

### Motorized panoramic imaging

Dimension + Multiposition

Standard + Multiposition

With cellSens Dimension or Standard software, wide area imaging using a motorized stage is fully automated with the optional Multiposition solution. When combined with Dimension and a motorized Z, this function can correct for the effects of sample distortion and tilting.




P1

P2

P3

P4




Panoramic imaging



## Well Plate Navigator

Dimension + Multiposition + Well Plate Navigator

The Well Plate Navigator automatically scans and acquires images from standard and customized well plate formats. All acquired images, sample positions, and user comments can be saved into a structured database for rapid centralized access. Move to the center of any well in a single click. Wells can be selected individually, by row or column, or in discontinuous groups. Apply unique multidimensional acquisition settings to a single well or multiple selected wells in one step. The Well Plate Navigator can execute multiple experiments within a single well plate in support of complex experiments.



## Extended Focus Imaging

Dimension  
Standard + Manual Process

Create a single in-focus image from successive image planes as the focus knob is turned using the Extended Focus Imaging (EFI) function. A motorized focus drive fully automates EFI acquisition. EFI composite images can also be created directly from previously captured Z-stacks.



Original image


Extend Focus Image

## Viewing and Processing

### Deconvolution

Dimension + CI Deconvolution

Olympus' optional constrained iterative (CI) Deconvolution Solution employs robust CI MLE algorithms to produce improved resolution, contrast, and dynamic range, with industry-leading speed. Each cellSens Dimension software license includes the most widely requested deblurring techniques such as 2D deconvolution, nearest neighbor, and Wiener filter.



Original image

Deconvolved image

Cell line: Human cervical cancer cell line HeLa cell  
Immunostaining: Hec1 staining (green, Alexa Fluor 488),  $\alpha$ -tubulin staining (red, Alexa Fluor 568), DAPI staining (blue)  
Mitotic HeLa cell derived from human cervical cancer.  
Mitotic spindle and kinetochores are stained with anti- $\alpha$ -tubulin (red) and anti-Hec1 (green) antibodies, respectively. Chromosomes interact with microtubules constituting mitotic spindle via kinetochores, protein structure assembled on centromere region of chromosomes.

#### Image data courtesy of:

Department of Molecular Oncology, Institute of Development, Aging and Cancer, Tohoku University  
Masanori Ikeda and Kozo Tanaka

### Solutions

Each cellSens software package can be expanded for specific applications by adding on the available "Solutions" noted below

Dimension		
CI Deconvolution	Multiposition	Well Plate Navigator
Count & Measure	Ratio/FRET	Database Core
Database Client	NetCam	Photo Manipulation
Standard		
Multichannel Acquisition	Multiposition	Count & Measure
Manual Process	Database Core	Database Client
NetCam	Confluency Checker	
Entry		
Database Client		


# A Range of Easy-to-Use Functions Turn Your Findings into Compelling Presentations

## Measurement and Analysis

### Manual Measurement

- Dimension
- Standard
- Entry


Distances between points, areas, intensity measurements, and morphological parameters are accessible using the cellSens software measurement tools. Measurement data is saved as an image layer that can be exported to MS Excel (except cellSens Entry) and cellSens workbook formats, or viewed using OlyVIA, a free image viewer software package.



### Automatic Object Measurement and Classification

- Dimension + Count & Measure
- Standard + Count & Measure

Count & Measure adds object detection for automated nuclei counting and classification. This solution extends the existing set of manual measurements in the cellSens software. Perform automatic object measurement and classification easily, using an interactive object based analysis that automatically links objects to their measurements.




### Intensity Analysis


- Dimension

Graphically depict intensity and ratio values defined by Regions of Interest (ROIs) and adjust ROI placement to compensate for cell movement. Export data directly to MS Excel. Convert variations of intensity to hue and brightness using Intensity Modulated Display (IMD) to visually enhance the fine image structure within ratio or FRET. The Ratio/FRET solution is used to display and analyze real-time ratiometric imaging and data. FRET analysis of both sensitized emission and acceptor photo-bleaching is also supported within this user friendly workflow.

The Photo-Manipulation solution can be used for the curve-fitting analysis of FRAP images.



Ratio image display/analysis




FRAP analysis function

### Object tracking



- Dimension + Count & Measure + Tracking

In time lapse imaging, moving objects can be automatically detected, tracked, and analyzed. cellSens Tracking provides a powerful and intuitive tool to quantify dynamic processes such as cell movement and division.



Time-dependent Change in Intensity of Cells

Image data courtesy of:  
Kazuhiro Yagita, M.D. Ph.D. Department of Physiology and Systems Bioscience,  
Kyoto Prefectural University of Medicine


## Collaboration and Communication

### Reporting

Dimension

Easily drag-and-drop image property data, measurement data, and user-customized fields into a report template using the convenience of a built-in tool to produce reports in MS Word.\* Collaborate with colleagues and communicate results quickly and easily.

\*Requires Word version 2007 or later



### Remote Live Image


Dimension

+ NetCam

Standard

+ NetCam

The cellSens NetCam solution facilitates the transfer of live or static imaging over a network for teaching, mentoring, or supervision.



### Database

Dimension

+ Database Core

or Database Client

Standard


+ Database Core

or Database Client

Entry

+ Database Client

The Database Core solution enables users to create shared, user-definable databases with full control over user access. The database stores images, associated image properties, user comments, and any other related files that a user wishes to include. The interactive query tool makes it easy to find the data and provides automatic previews of each queried image. Conveniently read and write to a shared database from different stations with the Database Client solution.



### Combination of Database and Well Plate Navigator

Dimension


+ Multiposition

+ Well Plate Navigator

+ Database Core

or Database Client

In combination with the Well Plate Navigator solution, the Database solution greatly improves the efficiency of viewing and analyzing well plate images with a large amount of data. By clicking on icons for image information such as the date, file name, or well plate number, any selection of captured images can be viewed for further analysis. This solution also enables users to view captured images and continuously analyze selected images (the Batch Macro function) via the well plate GUI.





## cellSens functions

		✓ : included functions / : optional solutions		
		DIMENSION	STANDARD	ENTRY
Layout	User experience customization	✓		✓
	Overlay multiple images	✓	✓	
View	Document groups for side-by-side image comparison	✓		✓
	Movie playback	✓	✓	✓
	Tile view (multiple images in a single data set shown side-by-side)	✓	✓	✓
	Slice view for orthogonal plane viewing of 3D or time-lapse data sets	✓		
	Voxel viewer for isosurface and volumetric rendering of 3D and 4D data sets	✓		
Image Acquisition	Snap/movie acquisition	✓	✓	✓
	Time-lapse at specified interval	✓		
	Automated multi-wavelength	✓	Multi-channel Acquisition	
	Z-stack	✓		
	Multi-dimensional (XYZT and wavelength)	✓		
	Graphical Experiment Manager	✓		
	Manual panoramic imaging (Instant MIA and Manual MIA)	✓	Manual process	
	Multiposition visitation and stage navigator	Multiposition	Multiposition	
	Automated panoramic imaging (auto MIA, requires motorized stage)	Multiposition	Multiposition	
	Instantly create EFI image (manual or motorized Z)	✓	Manual process	
	Simultaneous Multi-Color Imaging (requires two identical cameras or image splitter)	Ratio/FRET High-End Device	or	
	Live deblurring	✓		
High Dynamic Range Imaging (HDRI)	✓			
Multi-Well Plate Acquisition	Multiposition Well Plate Navigator	and		
Image Processing	Geometry/combine/filter processing	✓	✓	
	Fluorescence/Brightfield unmixing	✓		
	Deblurring (No/Nearest Neighbor, Wiener Filter)	✓		
	Kymograph	✓		
	2D deconvolution	✓		
Image Analysis	3D deconvolution (constrained iterative deconvolution)	CI Deconvolution		
	Region and line measurements	✓	✓	
	Phase analysis	✓		
	Object analysis and classification	Count & Measure	Count & Measure	
	Interactive measurement	✓	✓	✓
	Intensity plot over time/z	✓		
	Colocalization	✓		
	Object counting (Manual)	✓	✓	
	Object tracking	Count & Measure Tracking	and	
	Online ratio and kinetics	Ratio/FRET		
	Ratio analysis (off-line)	✓		
	FRET analysis	Ratio/FRET Life science analysis	or	
FRAP analysis	Photo manipulation Life science analysis	or		
Documentation and Collaboration	Cell count and confluency measurements	✓	Confluency Checker	
	Automatically compose MS Word reports	✓		
	Database image and data management solution for microscopy	Database Core Database Client	Database Core Database Client	Database Client
Remoting	Open database and load records/documents from database	Database Core Database Client NetCam	Database Core Database Client NetCam	Database Client
	Remote live image viewing			

## Products with confirmed functionality

		* 3 points angle, 4 points angle, arbitrary line, closed polygon, polyline, and perpendicular line only.			
		DIMENSION	STANDARD	ENTRY	
Olympus	Camera	DP21, DP22, DP26, DP27, DP72 <sup>1</sup> , DP73 <sup>2</sup> , DP74 <sup>4</sup> , DP80 <sup>2</sup>	✓	✓	✓
	Microscope	BX43, BX53, BX63, BX61, BX61WI, IX83, IX73, IX81, SZX16A	✓	✓	
	Peripherals	IX81-ZDC, IX81-ZDC2, IX3-ZDC, IX3-ZDC2	✓		
	Motorized XY stage	BX-DSU, IX3-DSU, IX2-DSU, U-CBF		Multiposition	
OSIS	Camera	XM10, XC10, XC30, XC50, UC30, UC50, UC90 <sup>3</sup> , LC20, LC30, SC30, SC50, SC100, SC180	✓	✓	✓
	Peripherals	cel*TIFF (multi-line, single line), MT20, USB-ODB converter, Real Time Controller (U-RTC and U-RTCE), U-FCB, U-STC, IX3-FRAP	✓		

## 3rd party products with confirmed functionality

				Note
Andor		iXon series, Zyla series, Neo		Requires High-End Camera solution
Hamamatsu		ORCA series, ImageM series, C11440-36U		
Photometrics		CoolSNAP HQ2, Evolve 512 Delta, Prime, Prime 95B		Several cameras require High-End Camera solution
Qimaging		MicroPublisher 3.3 RTV/5 RTV, Exi series, Qi Click, Retiga series, OptiMOS, Rolera Thunder		
Jenoptik		ProgRes C3, ProgRes C5		Available in cellSens Standard and Dimension
Prior		ProScan I, II, III, Optiscan II, III		
Ludl		MAC6000		
Objective Imaging		Oasis 4i		Available in cellSens Dimension, and requires Multi Position Solution for motorized stage use
Märzhäuser		Tango, Pilot stage		
Applied Scientific Instrumentation		MS-2000		
Vincent Associates		Uniblitz shutter (VCM-D1, VMM-D1, VMM-D3)		Available in cellSens Standard and Dimension
CoolLED		pE-1, pE-2, pE4000		
Excelitas		X-Cite 120 PC, X-Cite exacte, X-Cite XLED1, X-Cite 110LED, X-Cite 120LED, X-Cite TURBO		
Lumencor		SOLA SEII, SEII 365, Spectra X		Available in cellSens Dimension
Sutter		Lambda 10-3/10-B, Lambda DG4		
National Instruments		NI USB-6501		
Yokogawa		CSU-X1		Requires High-End Device solution

## Compatible image formats

Read and write	JPEG, JPEG2000, TIFF, BMP, AVI, PNG, VSI
Read only	GIF, PSD (Adobe PhotoShop), TIFF (DP-BSW, FSX100, MetaMorph), OIF/OIB/OIR (FLUOVIEW format), Cell, STK (MetaMorph), MRC (Medical Research Council)

## Recommended system requirements

OS*	Microsoft Windows 10 Pro (32-bit/64-bit) Microsoft Windows 8.1 Pro (32-bit/64-bit) Microsoft Windows 7 Ultimate/Professional (32-bit/64-bit) with SP1
OS Language	English, Simplified Chinese, Japanese, German, Russian (Entry and Standard) and Italian (Entry and Standard)
RAM /HDD / DVD drive	4GB or more is recommended / 1GB for installation / DVD-R DL compatible

\*cellSens Dimension and Dimension Desktop are only compatible with 64-bit operating systems.  
\*See detailed information: <http://www.olympus-lifescience.com/en/software/cellsens/>

## Image data courtesy of:

Hiroo Ueno, Ph.D. Department of Stem Cell Pathology, Kansai Medical University  
(cover page)

- OLYMPUS CORPORATION is ISO14001 certified.
- OLYMPUS CORPORATION is ISO9001 certified.

- All company and product names are registered trademarks and/or trademarks of their respective owners.
- Images on the PC monitors are simulated.
- Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.

[www.olympus-lifescience.com](http://www.olympus-lifescience.com)

**OLYMPUS**<sup>®</sup>

OLYMPUS CORPORATION  
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan

Printed in Japan M1718E-102017