The X8CAM4K Series HDMI/NETWORK/USB Multi-outputs CMOS Camera Help Manual

Со	ntents		I		
1	X8CAM4K Series Camera Application				
2	X8CAM4K Series Camera Datasheet and Functions (2)				
3	Dimensio	n of X8CAM4K Series Camera	4		
4	X8CAM4	Series Camera Packing Information	5		
5	Software	and App	6		
6	X8CAM4	Series Camera Configurations	7		
	6.1 0	amera working standalone with built-in XCamView software	7		
	6.2 0	onnecting camera to computers with USB3.0 port	8		
	6.3 0	amera working in WiFi mode (AP mode)	9		
	6.4 C	onnecting camera to the PC with LAN port	11		
	6.5 C applicatio	onnecting multi-cameras to the router through the LAN port/ WiFi STA mode for the non	etwork 12		
7	Instructio	ons for Using X8CAM4K Series Cameras with Different Microscopes	16		
	7.1 X	8CAM4K22MPA camera for use with biomicroscope via M42 microscope adapter	16		
	7.2 X	8CAM4K16MPA camera with trinocular stereo microscope	17		
	7.3 X	8CAM4K16MPA camera with metallographic microscope	17		
	7.4 X	8CAM4K22MPA with Canon lenses via M42 to EF mount adapter	18		
8	Brief Intr	oduction of X8CAM4K UI and Its Functions	19		
	8.1 XCamView UI				
	8.2 1	he camera control panel on the left or right side of the video window	19		
	8.3 T	he Measurement Toolbar on top of the video window	20		
	8.4 l	cons and functions of the Synthesis Camera Control Toolbar at the bottom of the video windc	w21		
	8.4.	L Browse	21		
	8.4.	2 Settings>Network	22		
	8.4.	3 Settings>Measurement	25		
	8.4.	1 Settings>Magnification	26		
	8.4.	5 Settings>Image Format	26		
	8.4.	5 Settings>Video	26		
	8.4.	7 Settings>Storage	27		
	8.4.	3 Settings>Files	27		
	8.4.	9 Settings>Time	28		
	8.4.	10 Settings>ISP	28		
	8.4.	11 Settings>EDF	28		
	8.4.	12 Settings>Stitch	29		
	8.4.	13 Settings>Voice Control	29		
	8.4.	14 Settings>Language	30		

	8.4.15	Settings>Miscellaneous	30
9	Sample Photos	Captured with X8CAM4K Series Camera	32
10	ToupTek ^{®.}	Contact Information	38

1 X8CAM4K Series Camera Application



Figure 1 The X8CAM4K Series Camera

The X8CAM4K series cameras have high resolution, high frame rate, and excellent image quality, making them suitable for professional applications that require extremely high levels of detail capture and colour reproduction. The basic characteristic is listed as below:

- Sony back-illuminated large target surface CMOS sensor
- With 10-bit ISP processing, it offers better color reproduction, remarkable sharpening and 3D noise reduction effects, and more accurate ROI white balance
- 4K HDMI/ NETWORK/ USB multiple video synchronous outputs
- 4K/1080P auto switching according to monitor resolution
- SD card/USB flash drive for captured image and video storage, support local preview and playback
- Support the capture and display of RAW format images
- Support Image Auto Upload to the server over the network
- Supports USB voice control module, enabling real-time control of the camera through voice commands
- for taking photos, recording videos, freezing, and other operations
- Provide real-time video EDF function and real-time video WDR output function
- Provide real-time Stitch function to obtain higher quality images through real-time processing
- Provide two sets of default ISP parameters for biological microscope and stereo microscope
- New UI interface, the browsing interface provids rich file operation functions, image to image comparison, image to real-time video comparison, multi-image EDF function, multi-image Stitch function

• Embedded XCamView for the control of the camera and image processing, supporting automatic edge finding and measurement functions

- The text input box supports the input of both Chinese and English languages
- ToupView/ToupLite software for PC
- iOS/Android applications for smart phones or tablets

2	X8CAM4K Series Camera Datasheet and Functions (2)
---	---

Order Code	Sensor & Size(mm)	Mount	Pixel(µm)	G Sensitivity Dark Signal	Sensor Output (FPS/Resolution)	Binning	Exposure(us)
X8CAM4K22MPA	IMX571(C) 1.7"(23.46x13.21)	M42*0.75	3.76x3.76	485mv with 1/30s 0.07mv with 1/30s	37@6240*3512	1x1	7~10 ⁶
X8CAM4K16MPA	IMX383(C) 1/1.1"(13.06x7.34)	С	2.4x2.4	470mv with 1/30s 0.21mv with 1/30s	60@5440*3060	1x1	41~10 ⁶

Camera Model	Video Saving(FPS/Resolution)	Picture	H DMI2.0(FPS/Resolution)	USB3.0(FPS/Resolution)	NETWORK(FPS/Resolution)
X8CAM4K22MPA	37@3840*2160	6240*3512	30@3840*2160 30@1920*1080	15@6240*3512 37@2688*1512 37@1920*1080	30@3840*2160 37@1920*1080 37@1280*720
X8CAM4K16MPA	60@3840*2160	5440*3060	60@3840*2160 60@1920*1080	25@5440*3060 45@2688*1512 60@1920*1080	30@3840*2160 60@1920*1080 60@1280*720



Figure 2 Available Ports on the Back Panel of the Camera Body

Interface or Button	Function Description
USB Mouse	Connect USB mouse for easy operation with embedded XCamView software
USB3.0	Connect USB flash drive to save pictures and videos Connect SG WiFi module to transfer video wirelessly in real time Connect USB microphone to record audio and video Connect USB voice control for enable real-time control of camera snap, recording, freezing, and other operations
USB Video	Connect PC or other host device to realize video image transmission
HDMI	Comply with HDMI2.0 standard. 4K/1080P format video output and supporting automatic switch between 4K and 1080P format according to the connected monitors
LAN	LAN port to connect router and switch to transfer video
SD	SD card slot, comply with SDIO3.0 standard and SD card could be inserted for video and images saving
ON/OFF	Power switch
LED	LED status indicator
DC12V	Power adapter connection (12V/1A)
Video Output Interface	Function Description
HDMI Interface	Comply with HDMI2.0 standard; 30fps@4K or 30fps@1080P (X8CAM4K22MPA) 60fps@4K or 60fps@1080P(X8CAM4K16MPA)
LAN Interface	Support real time resolution switching(4K/1080P/720P) H264 encoded video DHCP configuration or manual configuration Unicast/multicast configuration
WiFi Interface	Connecting 5G WiFi adapter (USB3.0 slot) in AP/STA mode
USB Video Interface	Connecting USB Video port of PC for video transfer MJPEG format video
Other Function	Function Description
Video Saving	Video format: 8M(3840*2160) H264 encoded MP4 file Video saving frame rate: 37fps(X8CAM4K22MPA); 60fps(X8CAM4K16MPA);
Image Capture	22M (6240*3512, X8CAM4K22MPA) JPEG/TIFF/RAW image in SD card or USB flash drive 16M (5440*3060, X8CAM4K16MPA) JPEG/TIFF/RAW image in SD card or USB flash drive (Default SD card priority, priority can be modified in settings)
Measurement Saving	Measurement information saved in different layer with image content Measurement information is saved together with image content in burn in mode
ISP	Exposure(Automatic / Manual Exposure) / Gain, White Balance(Manual / Automatic / ROI Mode), Sharpening, 3D Denoise, Saturation Adjustment, Gamma Adjustment, Contrast Adjustment, Brightness Adjustment, Hue Adjustment, Color to Gray, 50HZ/60HZ Anti-flicker, Dark Enhance, DRC Function

Image Operation	Zoom In/Zoom Out(Up to 10X), Mirror/Flip, Freeze, EDF, Stitch, Grids, Overlay, PIP, Browser(including Picture Browsing, Video Playback, Video Compare, Picture Compare, EDF, Stitch, Image Processing), Measurement Function				
Embedded RTC(Optional)	To support accurate time on board				
Restore Factory Settings	Restore camera parameters to its factory status				
Multiple Language Support	English / Simplified Chinese / Traditional Chinese / Korean / Thailand / French / German / Spanish / Japanese / Italian / Russian / Dutch / Portuguese				
	Software Environment under Network/USB Video Output				
White Balance	Auto White Balance				
Color Technique	Ultra-Fine Color Engine				
Capture/Control SDK	Windows/Linux/macOS/Android Multiple Platform SDK(Native C/C++, C#/VB.NET, Python, Java, DirectShow, Twain, etc)				
Recording System	Still Picture or Movie				
Operating System	Microsoft® Windows® XP / Vista / 7 / 8 / 8.1 / 10 / 11(32 & 64 bit) OSx(Mac OS X) Linux				
	CPU: Equal to Intel Core2 2.8GHz or Higher				
	Memory: 4GB or More				
DC De minemente	USB interface: USB 2.0 interface or higher				
PC Requirements	Ethernet Port: RJ45 Ethernet Port				
	Display:19" or Larger				
	CD-ROM				
Operating Environment					
Operating Temperature (in Centidegree)	-10°~ 50°				
Storage Temperature (in Centidegree)	-20°~ 60°				
Operating Humidity	30~80%RH				
Storage Humidity	10~60%RH				
Power Supply	DC 12V/1A Adapter				

3 Dimension of X8CAM4K Series Camera





4 X8CAM4K Series Camera Packing Information



Figure 5 X8CAM4K Series Camera Packing Information

	Standard Packing List				
Α	Gift box : L:25.5cm W:17.0cm H:9.0cm (1pcs, 1.7Kg/ box)				
В	X8CAM4K Camera				
С	Power Adapter: Input: AC 100-240V 50Hz/60Hz, Output: DC 12V 1A American standard: Model: POWER-U-12V1A(MSA-C1000IC12.0-12W-US): UL/CE/FCC European standard: Model: POWER-E-12V1A(MSA-C10001C12.0-12W-DE): UL/CE/FCC EMI standard: FCC Part 15 Subpart B EMS standard: EN61000-4-2,3,4,5,6				
D	USB Mouse				
Е	HDMI2.0 Cable				
F	USB3.0 A male to A male gold-plated connectors cable /1.5m				
	Optional Accessory				
G	Voice Control Module				
Η	SD Card(16G or above; Speed: class 10)				
Ι	USB flash drive(USB3.0)				
J	Ethernet cable				
K	USB WiFi adapter (Shape will vary with different models)				
L	ToupTek 1.2/1.5X Microscope Adapter (For X8CAM4K22MPA), ToupTek engineer will help you to determine the right microscope adapter for your application				

5 Software and App

The software or the APP can be downloaded from the following link:

Windows: <u>https://www.touptekphotonics.com.cn/download/?dlID=0</u> macOS: <u>https://www.touptekphotonics.com.cn/download/?dlID=1</u> Linux: <u>https://www.touptekphotonics.com.cn/download/?dlID=2</u> Android: <u>https://www.touptekphotonics.com.cn/download/?dlID=3</u> iOS: <u>https://www.touptekphotonics.com.cn/download/?dlID=4</u>

6 X8CAM4K Series Camera Configurations

You can use the X8CAM4K series camera in 5 different ways. Each application requires different hardware environment.

6.1 Camera working standalone with built-in XCamView software

For this application, apart from the microscope, you only need an HDMI monitor, the supplied USB mouse, and the camera embedded XCamView software. A computer or a network connection is not required to operate the camera in this application. The steps to start the camera are listed as below:



Figure 6 X8CAM4K Series Camera with the HDMI Monitor

Connect the camera to a HDMI monitor using the HDMI cable;



Insert the supplied USB mouse to the camera's USB Mouse port;



Insert the supplied SD card/USB flash drive into the X8CAM4K series camera SD card slot/USB3.0 slot;



Connect the camera to the power adapter and turn it on;



Turn on the monitor and view the video in the XCamView software. Move the mouse to the left, top or bottom of the XCamView UI, different control panel or toolbar will pop up and users could operate with the mouse at ease.



Figure 7 XCamView and X8CAM4K Series Camera in HDMI Mode

6.2 Connecting camera to computers with USB3.0 port

For Windows user (Windows XP (32bit), Windows 7/8/10/11 (32/64 bit)), please use ToupView.

For macOS and Linux user (macOS 10.10 or above or Linux distributions with kernel 2.6.27 or higher), please use ToupLite. The steps to start the camera are listed below:

Start the camera according to Sec. 6.1. After the camera is running, connect camera to computer with USB cable. Please use "USB Video" slot, The upper left corner of the HDMI graphics interface displays "USB3.0 Mode" or "USB2.0 Mode", indicating that a connection has been established with the PC.



Install ToupView/ToupLite on your PC or install ToupView App on the mobile device; Run the software ToupView/ToupLite, clicking the camera name in the Camera List group to start the live video as shown in Figure 8.



Figure 8 ToupView and X8CAM4K Series Camera in USB Mode

6.3 Camera working in WiFi mode (AP mode)

Please make sure your PC is WiFi enabled.



Figure 9 The PC or Mobile Device Connect to the Camera through WiFi For Windows user (Windows XP (32bit), Windows 7/8/10/10/11 (32/64 bit)), please use ToupView.

For macOS and Linux user (macOS 10.10 or above or Linux distributions with kernel 2.6.27 or higher), please use ToupLite. When connecting the camera with a mobile device, the free ToupView App is required. Just make sure that the mobile device uses iOS 11 or higher/Android 5.1 or higher operating systems.

The steps to start the camera are listed below:

Start the camera according to Sec. 6.1. After the camera is running, move the mouse to the bottom of the GUI and clicking the 💥 button on the Synthesis Camera Control Toolbar at the bottom of the video window, a small window called Settings will pop up as shown below. Click Network> WiFi property page and choose the AP in the WiFi Mode edit box(The factory default configuration is AP mode).

X	Settings	\times
Network	General LAN WiFi Tool	
Measurement	WiFi Mode: AP	
Magnification	Channel: 36	
Image Format	Password: 12345678	
Video		
Storage		
Files		
Time		
ISP		
EDF		
Stitch		
Voice Control		
Language		
Miscellaneous		
	Close	pply

Plug the USB WiFi adapter into the camera's USB3.0 port, the upper left corner of the HDMI graphics interface will display "AP mode";



Install ToupView/ToupLite on your PC or install ToupView App on the mobile device, connect the PC or mobile device to the camera's WiFi AP point; The network name (SSID) and the WiFi password (The default one is 12345678) can be found on the camera's Setting>Network> WiFi page in AP mode.



Start ToupView/ToupLite software or ToupView App and check the configuration. Normally, the active X8CAM4K series cameras will be automatically recognized. The live image of each camera is shown in Figure 10. For the display, the Camera List group is used in ToupView/ToupLite software, and the Camera Thumbnail is used in ToupView App.



Figure 10 ToupView and X8CAM4K Series Camera in WiFi AP Mode

6.4 Connecting camera to the PC with LAN port

This application uses the camera as the network camera. User must configure the IP of the camera and PC manually and ensure their IP addresses in the same net. The subnet mask and gateway of the camera and PC must be the same.



Figure 11 Connecting theX8CAM4K Series Camera with Ethernet Cable to the PC

Start the camera according to Sec. 6.1 after the camera is running, clicking 💥 button on the Synthesis Camera Control Toolbar at the bottom of the video window(See Figure 7), a small window called Settings will pop up as shown below on the left side, clicking LAN property page, uncheck the DHCP item. Input IP Address, Subnet Mask and Default Gateway for the camera. Designate Internet Protocol Version 4 (TCP/IPv4) Settings page's IP address on the PC with similar configuration as shown below on the right side but with different IP address.

\$Y	Settings	Internet 协议版本 4 (TCP/IPv4) Properties ×
Network Measurement Magnification	General LAN WiFi Tool DHCP Unicast O Multicast IP Address: 192 168 100 2 2 100 2 100 2 100	General You can get IP settings assigned automatically if your network supports this capability. Otherwise, you need to ask your network administrator
Image Format Video Storage Files Time ISP	Subnet Mask: 255 255 0 Default Gateway: 192 168 100 1 MAC Address: D0:5F:64:10:00:00	O Obtain an IP address automatically Use the following IP address: IP address: 192 . 168 . 100 . 3 Subnet mask: 255 . 255 . 255 . 0
EDF Stitch Voice Control Language Miscellaneous		Default gateway: 192 . 188 . 100 . 1 Obtain DNS server addresse automatically Use the following DNS server addresses: Preferred DNS server: Alternate DNS server:
	Close Apply	Validate settings upon exit Advanced OK Cancel

Figure 12 Configure theX8CAM4K Series Camera IP

Figure 13 Configure the PC's IP

After the above configurations are finished, user can connect the X8CAM4K series camera to the computer through the Ethernet cable as shown below:

Connect the LAN port with the Ethernet cable to the PC's network port, the upper left corner of the HDMI graphics interface will display IP address;



Insert the supplied SD card/USB flash drive into the X8CAM4K series camera's SD card slot/USB3.0 slot;



Install ToupView/ToupLite on your PC or install ToupView App on the mobile device; Run the software ToupView/ToupLite, clicking the camera name in the camera list starts the live video as shown in Figure 10.

6.5 Connecting multi-cameras to the router through the LAN port/WiFi STA mode for the network application

In LAN/ WiFi STA mode, the camera connects to the router by LAN port/ WiFi STA mode. If a router with LAN/ WiFi capability is used, users could connect the router with Ethernet cable/ WiFi to control the camera.



Figure 14 Multi X8CAM4K Series Cameras Connecting to the Router through the LAN Port/ WiFi Style

The connection and configuration are just the same as in Sec.6.1 or Sec. 6.4. But here, users need to check DHCP. If Multicast is disabled or is not supported, users should only select Unicast. If Multicast is supported by the network, users could select Multicast to achieve a better performance, especially in the case that multi-users connecting to the same camera. In addition, please guarantee that the broadcasting function is enabled in the network.

Active X8CAM4K series camera is recognized by ToupView/ToupLite software or ToupView App and they are displayed as a camera list or thumbnail in the software or app as shown in Figure 10.

S.	Settings	×
Network Measurement Magnification Image Format Video Storage	General LAN WiFi Tool Ø DHCP ● Unicast ○ Multicast IP Address: ● Subnet Mask: ● Default Gatess: ● MACC Address: ●	×
Files Time ISP EDF Stitch Voice Control Language Miscellaneous		
	Close Apply	у

Or start the camera according to Sec. 6.1. After the camera is running, move the mouse to the bottom of the video window and clicking the window and clicking the button on the Synthesis Camera Control Toolbar at the bottom of the video window, a small window called Settings will pop up as shown below. Clicking Network> WiFi property page and choosing the STA in the WiFi Mode edit box(The factory default configuration is AP mode). Choice or input the to be connected router's SSID and Password as shown below:

Ľ	Settings	×
Network Measurement Magnification Image Format Video Storage Files Time ISP EDF Stitch Voice Control Language Miscellaneous	General LAN WiFi Tool WiFi Mode: STA ▼ SSID: ▼ □ Password: ▼	

Install ToupView /ToupLite software on your PC. Alternatively, install the free ToupView App on the mobile device;

Plug the Ethernet cable into the camera's LAN port and the other end to the PC (for those connected to router with LAN Port), the upper left corner of the HDMI graphics interface will display IP address;



Or plug the USB WiFi adapter into the camera's USB3.0 port(for those connected to router with WiFi STA mode), the upper left corner of the HDMI graphics interface will display "STA Mode";



Finally, as shown below, 2 X8CAM4K series cameras are connected to the router with LAN cable and 2 X8CAM4K series cameras are connected to the same router with WiFi STA mode (The number of the cameras, the connection mode (LAN or WiFi STA) connected to the router are determined by the router performance).



Make sure that your PC or your mobile device is connected to the LAN or WiFi of the router; Start ToupView/ToupLite software or ToupView App and check the configuration. Normally, active X8CAM4K series cameras are automatically recognized. The live image of each camera is displayed. For the display, Camera List group is used in ToupView/ToupLite software, and Camera Thumbnail is used in ToupView App; Select theX8CAM4K series camera you are interested in. To do so, double click the camera's name in Camera List tool window if you use ToupView /ToupLite software; If you use ToupView App, tap the camera's thumbnail in Camera List page(See Figure 15)

About the routers/switches

It is suggested that routers/switches supporting WiFi 5G should be selected to achieve better wireless connection experience.



Figure 15 ToupView and X8CAM4K Series Camera in LAN port/WiFi STA mode

- 7 Instructions for Using X8CAM4K Series Cameras with Different Microscopes
- 7.1 X8CAM4K22MPA camera for use with biomicroscope via M42 microscope adapter



Figure 16 Biomicroscope + M42 Micro Adaptor +X8CAM4K22MPA



Figure 17 X8CAM4K22MPA Camera for Use with Biomicroscope Via M42 Microscope Adapter

7.2 X8CAM4K16MPA camera with trinocular stereo microscope



Figure 18 X8CAM4K16MPA Camera with Trinocular Stereo Microscope

7.3 X8CAM4K16MPA camera with metallographic microscope



Figure 19 X8CAM4K16MPA Camera with Metallographic Microscope

7.4 X8CAM4K22MPA with Canon lenses via M42 to EF mount adapter



Figure 20 X8CAM4K22MPA with Canon lenses via M42 to EF mount adapter



Figure 21 X8CAM4K22MPA with Canon Lenses Via M42 to EF Mount Adapter

8 Brief Introduction of X8CAM4K UI and Its Functions

8.1 XCamView UI

The X8CAM4K UI shown in Figure 7 includes a Camera Control Panel on the left of the video window, a Measurement Toolbar on the top of the video window and a Synthesis Camera Control Toolbar on the bottom of the video window.

	Notes				
1	To show the Camera Control Panel, move your mouse to the left or right of the video window. See Sec.8.2 for details				
	Move the mouse cursor to the top of the video window, a Measurement Toolbar will pop up for calibration and measurement operations. When				
2	user left-clicks the Float/Fixed button 📌 on the Measurement Toolbar, the Measurement Toolbar will be fixed. In this case the Camera				
	Control Panel will not pop up automatically even if users move mouse cursor to left or right side of the video window. Only when user left-				
	clicks the X button on the Measurement Toolbar to exit from measuring procedure will they be able to do other operations on the Camera				
	Control Panel, or the Synthesis Camera Control Toolbar. During the measuring process, when a specific measuring object is selected, an				
	Object Location & Attributes Control Bar $\land \forall \leqslant \geqslant \delta$ 🕅 will appear for changing location and properties of the selected object. See				
	Sec.8.3 for details.				
	When users move mouse cursor to the bottom of the video window, the Synthesis Camera Control Toolbar will pop up automatically.				
3	\oplus \bigcirc \bigotimes \bigotimes \bigotimes \bigcirc \bigotimes \bigotimes \bigoplus \bigoplus \bigoplus \bigoplus \bigoplus \bigoplus \bigotimes \bigoplus \bigotimes \bigotimes \bigoplus \bigotimes \bigotimes \bigoplus \bigotimes \bigotimes \bigoplus \bigotimes				

8.2 The camera control panel on the left or right side of the video window

The Camera Control Panel controls the camera to achieve the best video or image quality according to the specific applications; It will pop up automatically when the mouse cursor is moved to the left or right side of the video window (in measurement status, the Camera Control Panel will not pop up. The Camera Control Panel will only pop up when the measurement process is finished or terminated while user's cursor on the left edge of the video window). Left-clicking sutton to achieve Display/Auto Hide switch of the Camera Control Panel.

Camera Control Panel	Function	Function Description
	Snap	Capture image and save it to the SD card or USB flash drive
	Record	Record video and save it to the SD card or USB flash drive
	Auto Exposure	When Auto Exposure is checked, the system will automatically adjust exposure time and gain according to the value of exposure compensation
	Exposure Compensation	Available when Auto Exposure is checked. Slide to left or right to adjust Exposure Compensation according to the current video brightness to achieve proper brightness value
🗲 Camera Control Panel	Exposure Time	Available when Auto Exposure is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
	Gain	Adjust Gain to reduce or increase brightness of video. The Noise will be reduced or increased accordingly
Auto Exposure	Red	Slide to left or right to decrease or increase the proportion of Red in RGB on video
Exposure Compensation: 7	Green	Slide to left or right to decrease or increase the proportion of Green in RGB on video
Exposure Time: 16.667ms	Blue	Slide to left or right to decrease or increase the proportion of Blue in RGB on the video
Gain: 54	Auto	White Balance adjustment according to the window video every time the button is clicked
White Balance:	Manual	Adjust the Red, Green or Blue item to set the video White Balance
Red: 519	ROI	Check the ROI item will display a red ROI rectangle on the video window, drag it to the interested area will perform the White Balance according to the area video data
Blue: 444	One Push	Perform a global White Balance based on image conditions
One Push	Sharpness	Adjust Sharpness level of the video
Sharpness: 20	Denoise	Slide left or right to denoise the video
Denoise: 10	Saturation	Adjust Saturation level of the video
Saturation: 50	Gamma	Adjust Gamma level of the video. Slide to the right side to increase Gamma and to the left to decrease Gamma.
Contrast: 50	Contrast	Adjust Contrast level of the video. Slide to the right side to increase Contrast and to the left to decrease Contrast.
Brightness: 50	Brightness	Adjust Brightness level of the video. Slide to the right side to increase Brightness and to the left to decrease Brightness.
Hue: 50	Hue	Adjust Hue level of the video. Slide to the right side to increase Hue and to the left to decrease Hue.
Scence: Stereo V	DC	For DC illumination, there will be no fluctuation in light source so no need for compensating light flickering
	AC(50HZ)	Check AC(50HZ) to eliminate flickering caused by 50Hz illumination
	AC(60HZ)	Check AC(60HZ) to eliminate flickering caused by 60Hz illumination
	Scence	Select different default parameters according to the type of microscope
	Default	Restore all the settings in the Camera Control Panel to default values

8.3 The Measurement Toolbar on top of the video window

The Measurement Toolbar will pop up when moving mouse cursor to any place near the upper edge of the video window. Here is the introduction of the various functions on the Measurement Toolbar:

🗳 🛛 Visible Pixel	▼ NA	▼煮∡∧、	/>		>O\$0	$\odot \odot \odot$	൙ & ാ	⊡☆	Sum 7	产质	Х
	Figur	e 22 The Measu	ement To	olbar on t	he Unne	r Side of	the Vid	eo Wii	ndow		

Icon	Function
¢.	Float/ Fix switch of the Measurement Toolbar
✓ Visible	Show / Hide Measurement Objects
Pixel ·	Select the desired Measurement Unit
NA	Select Magnification for Measurement after Calibration
*	Object Select
4	Angle
\land	4 Points Angle
٠	Point (Point Counter)
	Arbitrary Line
\checkmark	3 Points Line
	Horizontal Line
	Vertical Line
\times	3 Points Vertical Line
11	Parallel
	Rectangle
\diamond	3 Points Rectangle
0	Ellipse
\odot	5 Points Ellipse
Θ	Circle
\odot	3 Points Circle
0	Annulus
0	3 Points Annulus
°	Two Circles and its Center Distance
a a a a a a a a a a a a a a a a a a a	3 Points Two Circles and its Center Distance
0	Arc
T	Text
$\hat{\Box}$	Polygon
S	Curve
um um	Scale Bar
\nearrow	Arrow
13	Execute Calibration to determine the corresponding relation between magnification and resolution, which will establish the corresponding relationship between measurement unit and the sensor pixel size. Calibration needs to be done with the help of a micrometer. For detailed steps of carrying out Calibration please refer to ToupView help manual.
AÛ	Auto Measurement: Two Points Parallel, Circle Detect, Annulus Detect, Rectangle Detect, Polygon
(\mathbf{F})	Export the Measurement information to CSV file(*.csv)
83	Measurement Setup
Û	Delete all the measurement objects
X	Exit from Measurement mode
A ♥ < > ▲	When the measurement ends, left-click on a single measuring object and the Object Location & Properties Control Bar will show up. User could move the object by dragging the object with the mouse. But more accurate movement could be done with the control bar. The icons on the control bar mean Move Left, Move Right, Move Up, Move Down, Color Adjustment and Delete.

1) When user left-clicks Display/Hide button \checkmark on Measurement Toolbar, Measurement Toolbar will be fixed. In this case Camera Control Panel will not pop up automatically even if moving the mouse cursor to the left edge of the video window. Only when user left-click the \times button on Measurement Toolbar to exit from the measurement mode will they be able to doing other operations on Camera Control Panel or Synthesis Camera Control Toolbar.

2) When a specific Measurement Object is selected during the measurement process, Object Location & Attributes Control Bar $\land \heartsuit \land \diamondsuit \land \diamondsuit \land \diamondsuit \land \diamondsuit \land \diamondsuit \land$ will appear for changing the object location and properties of the selected objects.

8.4 Icons and functions of the Synthesis Camera Control Toolbar at the bottom of the video window



Figure 23 The Synthesis Camera Control Toolbar on the Bottom of the Video Window

Icon	Function	Icon	Function
\oplus	Zoom In the Video Window	Q	Zoom Out the Video Window
\leq	Horizontal Flip	×.	Vertical Flip
G+G	Color/gray		Video Freeze
EDF	EDF	Stitch	Stitch
#	Grids	(Image Overlay
PIP	PIP		Browse images and videos in the SD Card
SS -	Settings	(j)	Check the Version of XCamView

The Browsing function, for detailed introduction, please refer to Section 8.4.1.

The 💥 Setting function, for detailed introduction, please refer to Sections 8.4.2 to 8.4.15.

8.4.1 Browse

Clicking the 🗇 to browse the dxf, images, videos, and other files saved on the SD card or USB flash drive, as shown in the following figure.



Figure 24 Browsing UI

There are two browsing modes: List mode and Thumb mode. The default is Thumb mode.

Right click on an empty area to create a new folder.

Double-click the search bar in the upper right corner of the mouse and enter the file name to perform the search. Right click on an image file to Copy, Cut, Rename, Delete, Video Compare, Picture Compare, EDF, Stitch and

Details. Clicking on a thumb to select the 1^{st} image, and clicking on another thumb to select the 2^{nd} image (or selecting 2 images with frame), then clicking the right mouse button to bring up the context menu and select Picture Compare to analyze and compare the two images (Four images can also be compared). Clicking on a thumb to select 3 (or box select 3) pictures focusing on different targets in the same scene, you can perform depth of field compositing on the selected pictures. Clicking on a thumb to select $2\sim32$ (or box select $2\sim32$) pictures, The selected images can be stitch in ascending order of the numerical numbers in the file name.

Right click on a video file to Copy, Cut, Rename, Delete, Video Compare, and view detailed information(Details).



Figure 25 Image Processing

Double-click the thumbnail of the picture with the left mouse button to open the picture, and then right-click the picture to Gray Scale, Invert, Binary, Saturation, Gamma, Brightness, Filter Color, Extract Color, Auto Level, Auto Contrast, Histogram, Histogram Equalization, Flip, and other image processing functions, and then after the processing is completed, you can choose reset to revert back to the original picture, and also you can choose save or save as in the lower sidebar of the picture. The description of each function is as follows:

Gray Scale	Choose Gray Scale command to convert a color image to a Gray Scale image
Invert	Choose Invert command to reverse the pixel values of the active image
Binary	Binary is a kind of gray level process. If the gray of the pixel is greater than the given threshold, the pixel's color will be changed into white. Otherwise, the pixel's color will be changed into black
Saturation	Adjust the Saturation of the image
Gamma	Adjust the Gamma of the image
Brightness	Adjust the Brightness of the image
Filter Color	Choose Filter Color command to filter a special color channel from a color image. Select either Red, or Green, or Blue color to filter. For every pixel, if select Red color to filter, only information about the Red channel will be discarded, the Green and Blue information will remain there.
Extract Color	Choose Extract Color command to extract a special color channel from a color image. Select either Red or Green, or Blue color to extract. For every pixel, if selecting Red color to extract, only information about the Red channel will be kept, the Green and Blue information will be discarded.
Auto Level	The Auto Level command moves the level's sliders automatically to set highlight and shadow. It defines the lightest and darkest pixels in each color channel as white and black and then redistributes the pixels' color values proportionately
Auto Contrast	The Auto Contrast command automatically adjusts the overall contrastin an RGB image
Histogram	Used to show the distribution of brightness, R, G, B of an image over an image
Histogram Equalization	Used to improved image contrast
Flip	Flip image Horizontally/Vertically

8.4.2 Settings>Network

8.4.2.1 Settings>Network>General



Figure 26 Comprehensive Network General Settings Page

Name	The current camera name recognized as the network name

8.4.2.2 Settings>Network>LAN



Figure 27 Comprehensive Network LAN Settings Page

DHCP	Dynamic host control protocol allows DHCP server to automatically assign IP information to the camera. Only in Sec 6.4 LAN networking this item should be checked, so that cameras can automatically get IP information from routers/switches to facilitate networking operation;
Unicast/Multicast	By default, unicast function is used. Only in Sec 6.4 networking environment, when the router/switch has multicast function, camera can switch to multicast mode, which can save the network bandwidth consumed by the camera and facilitate the connection of more cameras in the same network;
IP Address	Every machine on a network has a unique identifier. Just as you would address a letter to send in the mail, computers use the unique identifier to send data to specific computers on a network. Most networks today, including all computers on the Internet, use the TCP/IP protocol as the standard for how to communicate on the network. In the TCP/IP protocol, the unique identifier for a computer is called IP address. There are two standards for IP address: IP Version 4 (IPv4) and IP Version 6 (IPv6). All computers with IP addresses have an IPv4 address, and many are starting to use the new IPv6 address system as well. Users must manually configure their IP addresses on the camera side and computer side. The IP addresses set on the camera side and computer side should be in the same network segment. The specific settings are shown Figure 28. It's usually a private address. Private address is a non-registered address used exclusively within an organization. The internal private addresses retained are listed below: Class A 10.00-10.255.255; Class B 172.16.0-172.31.255.255; Class C 192.168.0-192.168.255.255. The suggested IP address is Class C.
Subnet Mask	Subnet Mask is used to distinguish network domain from host domain in 32-bit IP address;
Default Gateway	A default gateway allows computers on a network to communicate with computers on another network. Without it, the network is isolated from the outside. Basically, computers send data that is bound for other networks (one that does not belong to its local IP range) through the default gateway; Network administrators configure the computer's routing capability with an IP range's starting address as the default gateway and point all clients to that IP address.
MAC Address	Camera-independent physical address that identifies the network device.

Uncheck the DHCP and select the Unicast item, user still need to set the IP address, Subnet mask and Default Gateway as shown below:

X	Settings	×
Network	General LAN WiFi Tool	1
Network Measurement Magnification Image Format Video Storage Files Time EDF Stitch Voice Control Language Miscellaneous	General LAN WiFi Tool □ DHCP • Unicast ○ Multicast IP Address: 192 168 100 2 Subnet Mask: 255 255 0 0 Default Gateway: 192 168 100 1 MAC Address: D0:5F:64:10:00:00 1	
	Close Appl	y

Figure 28 Manual DHCP and Unicast

Uncheck the DHCP and select the Multicast item, user still need to set the IP address, Subnet Mask and Default Gateway as shown below:

\$	Settings	×
Network	General LAN WiFi Tool	
Measurement Magnification mage Format Storage Files Files Files SP EDF Stitch Voice Control .anguage Miscellaneous	Current Law Will Tool □ DHCP ○ Unicast ● Multicast IP Address: 192 168 100 2 Subnet Mask: 255 255 255 0 Default Gateway; 192 168 100 1 MAC Address: D0:5F:64:10:00:00 1 MAC Address: D0:5F:64:10:00:00 1	
	C	ose Apply

Figure 29 Manual DHCP and Multicast

8.4.2.3 Settings>Network> WiFi

X			Settings		×
Network Measurement	General LAN WiFi Mode:	WiFi Tool		V	
Magnification	SSID:			V	
Image Format	D Password:				
Video					
Storage					
Files					
Time					
ISP					
EDF					
Stitch					
Voice Control					
Language					
Miscellaneous					
					Class. Analys

Figure 30 Network Setup

Wi-Fi Mode	AP/STA mode to select;
Channel/SSID	Channel for the AP mode and SSID for the STA mode. Choice or input the to be connected router's SSID. Here, the SSID is the router's SSID;
Password	Camera Password for the AP mode. Router Password for the STA mode

8.4.2.4 Settings>Network> Tool

X	Settings	×
Network	General LAN WiFi Tool	
Measurement Magnification Image Format Video Storage Files Time ISP EDF Stitch Voice Control Language Miscellaneous	Image Auto Upload Enable Server IP Address: Connect Connect Description The Image Auto Upload function requires the installation of Server software on the computer.	
	Close Ap	

Figure 31 Comprehensive Network Tool Settings Page

Image Auto Upload	Select whether to enable or not;	
	When the WiFi mode is in AP mode, ensure that the PC is connected to the camera's AP, open the Server, click Update, and the IP address assigned by the camera to the PC will be displayed. Ensure that the Server has enabled; Manually enter the IP address and port on the camera end and click Connect. The left corner of the interface will display "Connected to Server", indicating successful connection. Click the snap button/click the left mouse button/use an external device to snap, The Server will display the number of Detections and Total Downloads, indicating successful Image Auto Upload;	
Server IP Address	When the WiFi mode is in STA mode, ensure that both the PC and camera are connected to the router's WiFi; When connected via LAN, ensure that the PC and camera are on the same LAN, open the Server, click Update, and the IP address assigned by the camera to the PC will be displayed. Ensure that the Server has enabled; Manually enter the IP address and port on the camera end and click Connect. The left corner of the interface will display "Connected to Server", indicating successful connection. Click the snap button/click the left mouse button/use an external device to snap, The Server will display the number of Detections and Total Downloads, indicating successful Image Auto Upload;	
Port	Default 8888	
Connect	Ensure that the Server has enabled, click Connect, and the left corner of the interface will display "Connected to Server", indicating successful connection;	
Description: The Image Auto Upload function requires the installation of Server software on the computer.		
Note: Enable Image Auto Upload function, unable to save pictures to SD card or USB drive; If you need to save pictures to an SD card or USB drive, you need to first turn off the Image Auto Upload function. For detailed instructions on the Image Auto Upload function and the Server on the upper computer, please consult our company for more information		

8.4.3 Settings>Measurement

This page is used for the define of the Measurement Object properties.



Figure 32 The Measurement Setup

	Precision	Used for setting digits behind the decimal point for measurement results;
	Edge Detection	Select whether to enable the automatic edge search function and set the detection range;
Clabal	Auto Measurement	Used for define the level of accuracy used for auto measurement;
Global	Font Size	The font size of measurement data can be divided into three types: large, Middle, and Small;
	Cursor	Select whether the cursor is a cross or a single cross and set the color of the cross or single cross;
	Miscellaneous	Whether to hide the label when moving the measurement objects;
Calibration	Line Width	Used for defining width of the lines for calibration;

	Color	Used for defining color of the lines for calibration;
	EndPoint	Type: Used for defining shape of the endpoints of lines for calibration: Null means no EndPoint, rectangle means rectangle type of endpoints. It makes alignment more easily;
Point, Angle, L	Point, Angle, Line, Horizontal Line, Vertical Line, Rectangle, Circle, Ellipse, Annulus, Two Circles, Polygon, Curve	
	Left-click the 庄 along with the Measurement command mentioned above will unfold the corresponding attribute settings to set	
	the individual property of the Measurement Objects.	

8.4.4 Settings>Magnification

This page's items are formed by the Measurement Toolbar's Calibration command.



Figure 33 Comprehensive Magnification Settings Page

Name	Names such as 4X, 10X, 20X, 40X, 100X are based on magnification of the microscopes. For continuous zoom microscopes, ensure that the selected magnification coincides with the scale alignment line on the microscope zoom knob; Users could also edit the name of the magnification with other information, for example, microscope mode, users name, etc.
Resolution	Pixels per meter. Image device like microscopes have high Resolution value;
Clear All	Click the Clear All button will clear the calibrated magnifications;
Delete	Click Delete to delete the selected magnification;
Up	Select a row in the magnification and click Move Up to move up the currently selected magnification;
Down	Select a row in the magnification and click Move Down to move up the currently selected magnification;

8.4.5 Settings>Image Format

5	Settings
Setwork Measurement Magnification mage Format Storage Files Fime SP	Settings Image Format Image Format Image Format Image Format Image Format Image Format Image Format Image Format Image Format Image Format
DF titch /oice Control .anguage /liscellaneous	Layered Mode Layered Mode Measurement objects are saved in different layer with image data in the target file. User could edit the measurement objects in the target file with software on the PC.
	Close Apply

Figure 34 Comprehensive Image Format Settings Page

Image Format	JPEG: The extension of JPEG file can get very high compression rate and display very rich and vivid images by removing redundant images and color data. In other words, it can get better image quality with the least disk space. If measurement objects are available, the measurement objects will be burned into the image and the measurement cannot be edited. TIFF: TIFF is a flexible bitmap format mainly used to store images including photos and artistic images. RAW (Little-ending byte order,16bits, high bits filled with 0): RAW is an uncompressed and unprocessed image format that preserves all raw data directly obtained from the sensor of a digital camera.
Measurement Object Saving Method	Burn in Mode: The measurement objects are merged into the current image. User could not edit the measurement objects any more. This mode is not reversable. Layered Mode: The measurement objects are saved in different layer with current image data in the target file. User could edit the measurement objects in the target file with some software on the PC. This mode is reversable.

8.4.6 Settings>Video

X	Settings	
Network Measurement Magnification Image Format Video Storage Files Time EDF Stitch Voice Control Language Miscellaneous	Video Resolution 0 1280x720 0 1920x1080 © 3840x2160 Video Quality O Low O Middle © High Playback Fast Forward/Reverse Interval: 20 = seconds	
	Clos	se Apply

Figure 35 Comprehensive Setting of Video page

Video Resolution	Select a Video Resolution of 1280 x 720, 1920x1080 or 3840x2160;
Video Quality	Select Video Quality as low, medium, or high;
Video Playback	Fast Forward/Reverse internal in second unite for Video Playback

8.4.7 Settings>Storage



Figure 36 Comprehensive Setting of Storage Page

Preferred Storage	SD Card: Select it to save the video and image to the SD Card.
Page	USB Flash Drive: Select it to save the video and image to the USB Flash Drive.
	List the file system format of the current storage device
File System	FAT32: The file system of SD Card is FAT32. The maximum video file size of single file in FAT32 file system is 4G Bytes;
Format of the	exFAT: The file system of SD Card is exFAT. The maximum video file size of single file in FAT32 file system is 16E Bytes;
Storage Device	NTFS: The file system of SD Card is NTFS. The maximum video file size of single file is 2T Bytes.
	Unknown Status: SD Card not detected or the file system is not identified;
Note: For USB Flash Drive, USB 3.0 interface is preferred.	

8.4.8 Settings>Files

8		Settings	×
Network Measurement	Image File Name		
Magnification	O Auto	⊖ Manual	
Image Format	Prefix: IMG		
Video	√Video File Name		
Storage	🗆 Add Time Suffix		
	O Auto	🔿 Manual	
Time	Prefix: VID		
ISP			
EDF			
Stitch			
Voice Control			
Language			
Miscellaneous			
			Close Apply

Figure 37 Comprehensive Setting of Files Name

File	Name	
Paradigm		
Auto		With specified name as the Prefix and XCamView will add digital after the Prefix for the Image or Video file;
Manual		A file dialog will pop up to enter the Image or Video file name for the captured Image or Video.

8.4.9 Settings>Time



Figure 38 Time Setting

Time User can set Year, Month, Day, Hour, Minute and Second ital.in this page.	onth, Day, Hour, Minute and Second ital.in this page.

8.4.10 Settings>ISP



Figure 39 Comprehensive Setting of ISP Page

Auto Exposure	Define the maximum automatic exposure time;
Metering Mode	Select the Metering mode as the Central Weighted Average Metering, Evaluative Metering, Partial Metering, or Spot
Metering Mode	Metering;
WB ROI Color	Choosing the ROI rectangle line color and whether it is synchronized display as Camera Control Panel;
Clarity Factor	Select to display the clarity factor in the video window, otherwise the clarity factor will not be displayed;
Dark Enhance	Define the intensity value of Dark Enhance;
DRC	Define the intensity value of DRC;

8.4.11 Settings>EDF

*		Settings		\rightarrow
Network	Automatic Alignment-			
Measurement	O None	⊙ Shift	⊖ Shift + Scale	
Magnification	- Sencitivity			
Image Format	OHigh	Middle	OLow	
Video				
Storage	Window size	∩ Middla	● Small	
Files		Omidale	0 sman	
Time	Description			
ISP	Automatic alignment: S	olves offset issues of fused image	ges, but slows down the fusion process.	
EDF	image.	e detection accuracy of depth of	neid, perhaps reduce the quarty of fused	
Stitch				
Voice Control				
Language				
Miscellaneous				
			Close Ap	ply

Figure 40 Comprehensive Setting of EDF

Automatic Alignment	Optionally turn on auto-alignment when there is significant displacement or scaling between images;
Sensitivity	Select the sensitivity of EDF;
Window size	Select the window size for displaying real-time images during EDF;
Description	Automatic alignment: Solves offset issues of fused images, but slows down the fusion process. Sensitivity: Improves the detection accuracy of depth of field, perhaps reduce the quality of fused image.

8.4.12 Settings>Stitch



Figure 41 Comprehensive Setting of Stitch

Detection Precision	Define the level of detection precision;
Stitching Stride	Select the stitching stride;
Background Color	Select the background color of stitch;
Description	Detection Precision: The higher the detection precision, the more precision the stitching and the longer it takes Stitching Stride: The larger the stitching stride, the lower the precision and the faster the speed.

8.4.13 Settings>Voice Control

Ś	Setting	S	×
Vetwork Measurement Magnification Image Format Video Storage Files Files Files	Voice Control ☑ Enable Key Words: Snap freeze record/begin record	unfreeze end/end record	
SP EDF Stitch Voice Control			
			Close Apply

Figure 42 Comprehensive Setting of Voice Control

Key Words	Provide Key Words for "snap";
	Provide Key Words for "freeze";
	Provide Key Words for "record/begin record", "end/end record";
Note: After the camera i	s turned on, if the voice control module is not plugged in the Key Words information will not be displayed by default

8.4.14 Settings>Language

*	Settings	×
Network Measurement Magnification Image Format Video Storage Files Time ISP EDF Stitch Voice Control Language Miscellaneous	● English ○ Simplified Chinese (海体中文) ○ Traditional Chinese (海体中文) ○ Korcan (ピュマー) ○ Thailand (กษา Thu) ○ French (Francais) ○ German (Deutsch) ○ Spanish (Español) ○ Japanese (日本治) ○ Italian (italiano) ○ Russian (pycecuii) ○ Dutch (Nederlands) ○ Portuguese (Portugués)	
		aga Annly

Figure 43 Comprehensive Setting of Language Selection Setting Page

English	Set language of the whole software into English;
Simplified Chinese	Set language of the whole software into Simplified Chinese;
Traditional Chinese	Set language of the whole software into Traditional Chinese;
Korean:	Set language of the whole software into Korean;
Thailand	Set language of the whole software into Thailand;
French	Set language of the whole software into French;
German	Set language of the whole software into German;
Spanish	Set language of the whole software into Spanish;
Japanese	Set language of the whole software into Japanese;
Italian	Set language of the whole software into Italian;
Russian	Set language of the whole software into Russian;
Dutch	Set language of the whole software into Dutch;
Portuguese	Set language of the whole software into Portuguese;

8.4.15 Settings>Miscellaneous

	Settings
etwork Acasurement Angnification mage Format dideo Atorage illes ime SP SP DF titleh Sioie Control approace	Ruler Show Color: Ø Enable Overlay Support saving overlay information in Burn In Mode Grids Guids Support saving grids information in Burn In Mode Monitor Working Mode Show O
discellaneous	Carson Size: Middle Left Click Snap Right Click Record Camera Control Panel Display Location Left Right Click Record

Figure 44 Comprehensive Miscellaneous Settings Page

Ruler	Select to display the ruler in the video window, otherwise not to display the ruler. You can choose the ruler color;		
Measurement	Select to display the measurement toolbar in the video window, otherwise not to display the measurement toolbar;		
Overlay	Select to support saving graphics overlay information in fusion mode, otherwise it will not support;		
Grids	Select to support saving mesh information in fusion mode, otherwise not to support;		
Monitor Working Mode	Select to display the Monitor Working Mode in the video window, otherwise the Monitor Working Mode will not be displayed;		
Cursor	Choosing the Mouse size according to the screen resolution or personal preference; Select to Left Snap Right Record. If not selected, it will not Left Snap Right Record		
Camera Control Panel Display Location	Select the camera control panel to display on the left, right, or both sides of the HDMI interface;		

Camera	Parameters	Import the Camera Parameters from the SD Card or USB flash drive to use the previously exported Camera	
Import		Parameters	
Camera	Parameters	Export the Camera Parameters to the SD Card or USB flash drive to use the previously exported Camera Parameters	
Expon			
Reset to factory defaults		Restore camera parameters to its factory status;	

9 Sample Photos Captured with X8CAM4K Series Camera



Figure 45 Corn Stem.C.S Captured with X8CAM4K22MPA



Figure 46 Rhizopus.W.M Captured with X8CAM4K22MPA



Figure 47 Cross Section of Cotton Stigma Captured with X8CAM4K22MPA



Figure 48 FPC Captured with X8CAM4K16MPA



Figure 49 Carbide Tip Captured with X8CAM4K16MPA



Figure 50 Circuit Board Captured with X8CAM4K16MPA



Figure 51 Samples Photographed by the SOPTOP RX50M Metallographic Microscopy System Captured with X8CAM4K16MPA



Figure 52 Conductive Particles of the Liquid Crystal Screen Photographed by the DIC-100WVA Differential Interference Contrast (DIC) Microscopy System Captured with X8CAM4K16MPA

The X8CAM4K Series Camera Help Manual



Figure 53 Conductive Particles of the Liquid Crystal Screen Photographed by the DIC-100WVA Differential Interference Contrast (DIC) Microscopy System Captured with X8CAM4K16MPA



Figure 54 Metal Samples Photographed by the OUMIT CX43M Metallographic Microscopy System Captured with X8CAM4K16MPA



Figure 55 Semiconductor Samples Photographed by the OUMIT CX43M Metallographic Microscopy System Captured with X8CAM4K16MPA

10 ToupTek®--- Contact Information

	杭州图谱光电科技有限公司				
	杭州市西湖区西园五路6号奥强大厦1号楼15层				
	杭州, 310030, 浙江,				
•	中国				
	Hangzhou ToupTek Photonics Co., Ltd				
	15F, Aoqiang Building 1, No. 6, Xiyuan 5th Rd.,				
	Hangzhou, 310030, Zhejiang,				
	P.R.China				
	+86-571-8111-0735				
	+86-571-8111-0730				
S. 1	+86-571-8810-2638,				
	+86-18058780750 (手机/Mobile Phone)				
	FAX: +86-571-8668-3738				
	tphz@touptek.com				
\sim	Skype:	18058780750/ToupTek Photonics			
····)	QQ	2426878316			
	Wechat	18058780750			