


KoPa[®]



Innovative 5G WiFi (11ac) Interactive Microscope Classroom System

Install in minutes / Multiple operation systems supported / Free of tedious infrastructure

Support rapid digitization of pure optical trinocular microscope



Three interactive system solutions

1 One microscope vs. 1 - 8 smart devices. Suitable for small laboratory / class / workshop, etc.

| | | | |
|-------------------------|----------------|-----------------------|-----------------|
| 1080P FULL HD | 2.0MP 5.0MP | 4K ULTRA HD | 8.0MP 12.0MP |
|-------------------------|----------------|-----------------------|-----------------|




Solution:

- 5G WiFi (11ac) camera (choose from the following models) :
MC1210-W, MC810-W, MC510-W, MC210-W (please refer to the specifications of corresponding models)
- Windows PC software: Micro WiFi Lab
- Smart devices applications: Micro WiFi Lab
- Smart devices (purchased with a distributor or provide for oneself) : support iOS, Android and Windows systems;
- Hardware requires 5G WiFi (11ac), while laptop or desktop without 5G WiFi (11ac) requires an external network card.

2 One microscope vs. 11-50 smart devices. Suitable for large and medium-sized laboratories / classrooms / seminars and so on.

Solution:


- 5G WiFi (11ac) camera (choose from the following models) :
MC1210-W, MC810-W, MC510-W, MC210-W (please refer to the specifications of corresponding models)
- Master ViMatrix VM3500
- Windows PC software: Micro WiFi EDU
- Smart devices applications: Micro WiFi EDU
- Smart terminal (purchased with the agent or provide for oneself) : support iOS, Android and Windows systems;
- Hardware requires 5G WiFi 11ac, while laptop or desktop without 5G WiFi 11ac requires an external network card



3

12-60 microscope vs. 12-60 smart devices.

It is suitable for the microscope education activities of 12-60 people in universities, research institutes, hospitals and so on.



● Please refer to the following parts of this brochure for the system composition.

Example of 5G WiFi (11ac) Microscope Classrooms



LOWER COLUMBIA COLLEGE



Training




- No.1 Classroom : 49 sets
- No.2 Classroom : 25 sets
- No.3 Classroom : 25 sets
- No.4 Classroom : 17sets



Smartphones, tablets, and laptops are the carriers of modern networked knowledge. 5G WiFi (11ac) Interactive Microscope Classroom System gives us a great sharing platform that liberates you from device or location and give you unlimited access and function to knowledge -making learning fun and engaging us in the 21st technology. Let's face it, everybody that has a smartphone /tablet are walking computers with access to data on moments notice via text or verbal commands. The future is here. It only takes a microscope, a WiFi camera and an App to let student / teacher collaboration meet anywhere. Where the students take this technology from here is the most exciting part. We have successfully integrated every microscope in the classroom with this new technology.


System layout

Master 5G WiFi (11ac) ViMatrix





















Student-side

Real-time images under microscope are transmitted wirelessly to teacher-side computers and mobile mini PC via 5G WiFi (11ac) signals



System composition

- ✓ No tedious Infrastructure, efficient installation, plug and play.
- ✓ Individual smart devices are involved into education activities and the efficiency of study is improved exponentially.
- ✓ Not site restrictions, easy to move the system as you want.
- ✓ Support four major systems: Windows, Mac, iOS, and Android.
- ✓ Independent camera desinged for rapid upgrades in the future.
- ✓ Self-built LAN, high-speed and stable data transmission.

| | |
|---|---|
| <p>Optical microscope</p>   | <p>Teacher computer</p>  <p>Software: Micro WiFi EDU</p>   <p>When used on a PC, you need to install the interactive software Micro WiFi EDU which included in this System .</p> |
| <p>Optional 0.5x or 1x adapter</p>  | <p>5G WiFi (11ac) supported tablet and smartphone (please bring yours or purchase with distributors)</p>  <p>APP: Micro WiFi EDU</p>   <p>Scan provided QR code to download APP Micro WiFi EDU for iOS and android smart devices</p> |
| <p>Camera for microscope</p>   <p>Camera for student</p>  <p>Camera for teacher</p> | <p>Projection</p> <p>The teacher side computer can use the projector for large-screen display, but also can use the lecturing mode of this system, real-time stream such as microscopic images, comptuer documents to the smart devices of students, so as to obtain high-quality teaching effect.</p> |
| <p>Master 5G WiFi (11ac) ViMatrix VM3500</p>    | <p>Audio system</p> <p>Teachers and students can do one-to-one voice communication;to ensure audio quality, dedicated high-quality microphone for teachers is provided by designated agents;students can use their own headphones.</p> |
| <p>Student side Vimatrix Mini</p>    | <p>Classroom requirements</p> <p>Student desks requirements: please provide enough physical space to the student microscope and Vimatrix Mini, 3 student microscope cameras share one Vimatrix Mini. It is recommended that three student cameras have a physical distance of less than 2 meters . Two mains outlets should available (rated output voltage is AC 100-240V).</p> |



Olympus optical microscope



- Ergonomic grips add safety when retrieving the microscope from high places.
- Colored grips indicate the appropriate places to hold the microscope.
- Smooth, rounded design eliminates sharp edges.
- The CX23 microscope is one of the lightest in its class, with a total weight of approx. 5.9 kg (13.01 lbs).

1



Convenient and easily accessible power cable storage



Locking pin for easy binocular rotation

- A locking pin keeps the observation tube in place.
- Interpupillary adjustments, ranging from 48 to 75 mm, enables individual users to set it to their needs.
- Eyepoint adjustments accommodate the user's height and provide greater comfort.
- Left and right diopter adjustment enable optimal focus for each eye.

3




Interpupillary adjustment



Eyepoint height adjustment



Focus lock



Wide field of view

- The focus lock helps prevent the objectives from damaging the specimen.
- Field number (FN) of 20 provides a wide field of view.
- The CX23 microscope features plan achromat objectives that preserve outstanding image flatness throughout the field of view.
- The energy saving LED light source provides a long operation lifetime of 20,000 hours.
- Daylight balanced LED illumination preserves the vivid colors on HE-stained samples.

5



Angled arm enables comfortable carrying



Ergonomic grip for easy carrying

- A locking pin keeps the observation tube in place.
- Interpupillary adjustments, ranging from 48 to 75 mm, enables individual users to set it to their needs.
- Eyepoint adjustments accommodate the user's height and provide greater comfort.
- Left and right diopter adjustment enable optimal focus for each eye.

2



Fixed eyepieces



Rackless stage and stage cover



Inward-facing revolving nosepiece



Coaxial knobs

- An inward-facing rotating nosepiece facilitates a larger working area above the stage, enabling easier operation and specimen exchange or making it easier to apply oil to immersion objectives.
- Objectives with long working distances provide clearance for thick samples and slides.
- Bring specimens into focus quickly by using the coaxial knobs for fine or coarse adjustments. The microscope is built for durability and precise control, whether accessed from the left or right side.
- A torque release function helps prevent damage if excess force is applied to the coarse focusing knob or stage knob.

4



Cable storage compartment



Built-in security slot

- A storage compartment on the back of the CX23 microscope makes it easy to stow the power cable after use.
- Protect your investment with an optional custom designed wooden case.
- When the CX23 microscope needs to be left unattended, a built-in security slot enables attachment of an anti-theft cable.

6


Camera for microscope



Plug and play, powered by the ViMatrix Mini, safe and reliable. Do not change the mechanical structure of the microscope, do not destroy the original optical system of the microscope. Static and dynamic pixel are real synchronized, no pixel interpolation.

The camera only has a power button to prevent students from misoperation. Independent camera designed makes it easy to quickly upgrade a higher resolution camera in the future. When the camera needs to be repaired or replaced, it can be instantly unscrewed without waiting for the after-sales personnel to be present.

- ✓ 2.0MP
MC211-M
- ✓ 5.0MP
MC511-M
- ✓ 8.0MP
MC811-M
- ✓ 12.0MP
MC1211-M



Master 5G WiFi (11ac) ViMatrix VM3500

The unique ViMatrix controls all cameras, and the teacher-side PC has all administrator rights of the system.

- 9 strong signal antennas, can cover the diameter of up to 50 meters of spherical environment, to meet the requirements of teaching activities.
- Ultra-strong core to ensure stable real-time image and data transmission.
- Provides 1 WAN to support Internet connection and 4 LAN ports available to connect with teacher's PC.
- Multi-terminal data synchronization and transmission, the maximum transmission rate can reach 2600Mbps @ 5G Hz.




Student side 5G WiFi (11ac) ViMatrix Mini

- One high gain single frequency external antenna.
- Friendly-designed three 5V 2A USB charging interface, available for students to charge their smart devices in the classroom.
- Three power ports for microscope camera.
- Three cameras share one ViMatrix Mini.



Teacher Interactive software Micro WiFi EDU


Six main interactive teaching methods allows arbitrary composition of use



- User-friendly interface with icons, easy to identify and quickly use.
- Each function has been designed to be very accessible and helpful to education and research.
- Unique black background, eliminate miscellaneous light interference, improve microscope image recognition.

Scope of applications:

- | | | |
|---|---|---|
| <input checked="" type="checkbox"/> Education teaching | <input checked="" type="checkbox"/> Cell Lab | <input checked="" type="checkbox"/> Pathological Studies |
| <input checked="" type="checkbox"/> Blood Science | <input checked="" type="checkbox"/> Parasite Science | <input checked="" type="checkbox"/> Microbiology |
| <input checked="" type="checkbox"/> Material analysis | <input checked="" type="checkbox"/> Environmental studies | <input checked="" type="checkbox"/> Food/Beverage Control |
| <input checked="" type="checkbox"/> General cytology such as agriculture, forestry, livestock | | |



Lecturing mode



The teacher PC screen or the selected student smart device screen is streamed to all student smart devices synchronously.



Lecture materials such as PPT/ WORD/ EXCEL are displayed on the smart device of students in real time to enrich the teaching content.



The operation process of the teacher PC or the selected student smart device is displayed to all student smart devices synchronously.


Demonstrate student works

Live-demo the selected camera video to the entire student smart devices.



Monitor live-image under microscope camera

Teachers can easily capture all camera live-image, confirming that all cameras are connected successfully.




Annotate student works

Annotate and edit selected images and documents, and send the operation process to all student smart devices in real time.



Monitor smart devices screen

Real-time monitoring of all student smart devices screen to ensure that all students participate in learning. The students' annotations, measurements, and other operations in the APP are displayed on teacher's PC. The progress of the study can be checked at any time. The teacher can show the selected student smart devices screen to all students.



Multi-screen comparison teaching

- Static picture comparison
- Dynamic Video comparison
- Static vs. dynamic mixing comparison
- Support for cross-window annotating
- Support teacher-to-peer comparison with selected students

Teacher interactive software Micro WiFi EDU Interactive shortcuts

Multi-page camera monitoring window, capturing the teacher-side image and all student-side real-time image simultaneously.

Take picture, record, screen-capture, screen-record for single or multiple student-side images.

Smart settings for camera : restore all camera default parameters;
Restore all camera default names;
Disable parameters adjustment on all smart devices.
Modify the camera name; restore the current camera default name.

The teacher receives students' voice and answers one-on-one questions.

Delete all documents with one click and reset the system with one click (convenient for next class teachers and students).

Smart devices black and white list management: delete all white list and blacklist; allow Smart devices authorization message popup; The smart device is authorized to join into the system in order. Check all authorized smart devices list.

23 measurement tools to meet basic teaching requirements.




Teacher-side mobile Mini PC (optional)

The mobile Mini PC has the same control and functions as the teacher-side PC, enabling easy mobile teaching in the classroom.



Student smart device APP Micro WiFi EDU




Scan supplied "student seat QR code" to access micro image directly.

Take photos / videos and record important content in real time.


Annotate, measure, doodle, add pictures, improve learning efficiency.

The student smart devices and the teacher-side PC can send each other voice realizes individual counseling.

Generates learning / research reports within student APP easily and can sent to teacher's email immediately.



Tool bar



Annotation tools




System Image Samples


Taken by CX23 and camera MC811-M




Immunohistochemistry




Ovarian cancer membranes 40X



Colon pulp 40X




Appendix nucleus 40X




Cervical cancer nuclear slurry 40X

Insect




Wheel-worm 40X




Paramecium stons 40X


Pathological samples




40X



40X




40X




40X

Animals samples




Canine ileum 40X




Dog taste buds 40X

Plant



Pine young male cone longitudinal 10X



Lime tree stem cross 10X

Technical specifications

| | |
|--------------------------------|---|
| Microscope model | CX23 |
| Body | gnirevoc citsalp evitcetorp ,emarf latem tsac-eid munimulA |
| Optical system | UIS2 optical system (Infinity optical system) |
| Illumination system | Built-in transmitted illumination system, LED power consumption: 0.5 W (nominal value) |
| Focusing | Stage height movement (Coarse movement stroke:15mm) ,Coarse adjustment limit stopper <ul style="list-style-type: none"> • Torque adjustment for coarse adjustment knob • Fine focus knob (minimum adjustment gradations: 2.5 μm) |
| Revolving nosepiece | Fixed quadruple nosepiece with inward tilt |
| Stage | egats dexif lacinahcem tnevemov eriW Traveling range (X × Y): 76 mm × 30 mm, Specimen holder, Specimen position scale |
| Observation tube | Tube inclination 30° Interpupillary distance adjusting range: 48 - 75mm, Eyepoint adjustment: 370.0 - 432.9mm |
| Objectives | Plan Achromat, anti-fungal 4x NA: 0.10 W.D.: 27.8 mm 10x NA: 0.25 W.D.: 8.0 mm 40x NA: 0.65 W.D.: 0.6 mm 100xOil NA: 1.25 W.D.: 0.13 mm (CX23LEDRFS1 built-in) |
| Eyepiece(10x) | Field of view (FN): 20 (mildew-proof treatment) |
| Optional accessories | Reflector (CH20-mm), 15 X eyepiece (WHSZ15X-H: FN12, mildew-resistance treatment), Special wooden box, eyepiece micrometer, dark field diaphragm (CH2-DS+ CH2-FH) |
| Weight | Approximately 5.9 kg |
| Rated voltage/Electric current | AC100 – 240 V 50/60 Hz 0.4 A |
| Power consumption | Less than 2W |

| | MC211-M | MC511-M | MC811-M | MC1211-M |
|-----------------------|--|--|---|--|
| Student camera | MC211-M | MC511-M | MC811-M | MC1211-M |
| Teacher camera | MC211-W | MC511-W | MC811-W | MC1211-W |
| Resolution | 2.0MP | 5.0MP | 8.0MP | 12.0MP |
| Image sensor | SONY IMX291 CMOS | SONY IMX178 CMOS | SONY IMX274 CMOS | SONY IMX226 CMOS |
| Exposure mode | Rolling exposure | Rolling exposure | Rolling exposure | Rolling exposure |
| Maximum resolution | 1920 x 1080 (2,073,600 Resolution) | 2592 x 1944 (5,038,848 Resolution) | 3840 x 2160 (8,294,400 Resolution) | 4000x 3000 (12,000,000 Resolution) |
| Sensor size | 1/2.8" | 1/1.8" | 1/2.5" | 1/1.7" |
| Pixel size | 2.9 μm x 2.9 μm | 2.4 μm x 2.4 μm | 1.62 μm x 1.62 μm | 1.85 μm x 1.85 μm |
| Dynamic range | 128dB | > 80dB | > 80dB | > 80dB |
| Signal-to-noise ratio | 30dB | ≥ 50dB | ≥ 50dB | ≥ 50dB |
| Spectral response | 380-650nm | | | |
| Exposure capability | Real-time automatic, single-time automatic, manual adjustment | | | |
| White balance | Real-time automatic, single-time automatic, manual R B separately adjusted | | | |
| Record format | Photo Picture format: JPG Resolution: 1920x1080 | Photo Picture format: JPG Resolution: 2592x1944 | Photo Picture format: JPG Resolution: 3840x2160 | Photo Picture format: JPG Resolution: 4000x3000 |
| | Video Video format: MOV Resolution: 1920X1080@50FPS(50Hz), 1920X1080@60FPS(60Hz), 1280x720@50FPS(50Hz), 1280x720@60FPS(60Hz) | Video Video format: MOV Resolution: 2592X1944@30 FPS, 2560X1920@30 FPS, 2048x1536@30 FPS, 1920x1080@30 FPS | Video Video format: MOV Resolution: 3840X2160@25 FPS, 4096X2160@25 FPS, 2048x1536@25 FPS | Video Video format: MOV Resolution: 4000X3000@15 FPS, 4096X2160@25 FPS, 3840x2160@25 FPS, 2592x1944@25 FPS |


| Software operating requirements | |
|---------------------------------|---|
| System requirements | Microsoft® Windows® 7 / 8 /8.1/10(32 & 64 bit) |
| | CPU: i7 eight generation quad core or higher |
| | Memory: 8G or above |
| | Hard disk: at least 50 GB free space |
| | Graphics card: discrete graphics (NVIDIA display chip) |
| iOS requirements | Network card: 10/100/1000M adaptive |
| | iPhone X/ 8/ 8 Plus/7/ 7 Plus/6s/6s Plus/ 6/6 Plus iPad Pro(12.9 inches,2nd Generation)/(12.9 inches,1st Generation)/(10.5 inches)/(9.7 inches) iPad Air 2/iPad mini 4 and other new devices that support the 5G WiFi (11ac) protocol |
| Android requirements | iOS11.0 and later system |
| | Android 5.0 and later system |
| | CPU: Dual core 1.7GHz and later |
| | Memory RAM:3G or more |
| | Storage ROM: 8G or more recommended |
| | Support 5G WiFi (11ac) protocol |

| Master 5G WiFi (11ac) ViMatrix model | VM3500 |
|--------------------------------------|--|
| Interface | 4 x 10/100/1000BASE-T Ethernet RJ45 Interface 1 x 10/100/1000BASE-T Ethernet RJ45 Interface |
| Button | Power on/Off |
| Lights | Power Light, Status light |
| Antenna | 9 high-gain single-band antennas |
| External power supply | AC100~240V/1.5A (Max) |
| Wireless parameters | |
| Protocol supported | 5G WiFi IEEE802.11ac |
| Frequency | 5.180-5.825GHz |
| Signal rate | 2600Mbps |
| Transport security | WPA/WPA2, WPA-PSK/WPA-PSK2 encryption |
| WAN type | Dynamic IP |
| Protocol | Supports IPV4 and IPV6 |

| Student - side 5G WiFi (11ac) ViMatrix Mini model | WF01A |
|---|---|
| Interface | 3 x 10/100BASE -T Ethernet RJ45 Interface |
| Button | Power on/Off |
| Lights | Power Light, Status light |
| Antenna | 1 high-gain single-band antennas |
| External power supply | DC12V 5A |
| Wireless parameters | |
| Protocol supported | 5G WiFi IEEE802.11ac |
| Frequency | 5.180-5.825GHz |
| Signal rate | 450Mbps |
| Transport security | WPA/WPA2, WPA-PSK/WPA-PSK2 encryption |
| Protocol | Supports IPV4 and IPV6 |



Dimensions




About 5G WiFi (11ac)



5G WiFi (11ac) advantages over 2.4G WiFi

| 2.4G | 5G WiFi(11ac) |
|---|--|
| <p>Slow network speed, easy to be congested</p> <ol style="list-style-type: none"> 1. Less channels: The 11n protocol supports a maximum of 14 channels. 2. Channel bandwidth is narrow: The transmission bandwidth of each channel reaches to only 20M. 3. Ttransmission speed is slow: Through the existing technology is difficult to achieve high transmission speed. Current normal maximum transfer speed is only 300Mbps. | <p>Have a faster network and solve congestion</p> <ol style="list-style-type: none"> 1. More channels: The 11ac protocol supports more than 100 channels. 2. Larger channel bandwidth: The transmission bandwidth of each channel reaches up to 80M. 3. Higher transmission speed: MU-MIMO technology is very convenient to achieve a transmission speed of 2600Mbps, or even higher. |
| <p>Poor signal quality and heavy interference</p> <p>At present, the 2.4GHz band WiFi network is not only used on mobile phones, tablets, laptops, etc., but also various mobile devices such as handheld game consoles, bluetooth mouse/keyboard, bluetooth headsets, and microwave ovens. A large number of devices accumulated in a small frequency band and easily interfere with each other.</p> | <p>Better signal quality, less interference</p> <p>The 5G WiFi (11ac) band is less used and radio interference is greatly reduced, signal quality is greatly improved.</p> |

5G WiFi (11ac) supports more channel bands and has a faster network. The 5GWiFi (11ac)ViMatrix has a maximum transmission rate of 2,600 Mbps, easily covering a spherical range of 50 meters in diameter, and is suitable for teaching activities.



Difference between 5G communication and 5G WiFi (11ac)



5G communication refers to the fifth generation mobile communication technology, referred to as 5G, which is an extension after 4G, and its working frequency is 3.5GHz. 5G WiFi (11ac) refers to the fifth generation WiFi transmission technology, referred to as 5G WiFi. WiFi technology working in the 5.180-5.825 GHz band. The way they are implemented is completely different from the purpose. 5G communication is used for telephone calls, SMS communication, and mobile Internet access. The 5G WiFi (11ac) is used for large-scale and high-speed video and image transmission in small areas.

FC Verification No.: GZEM1801000017ITV



Verification No.: CSRT180084
Verification No.: CSTT180018

CE Verification No.: GZEM1801000016ITV



FCC IDENTIFIER: 2AFO3WF01A
CE Certification NO.: 18/0050/SZ



High-Tech Enterprise certificate number:
GR201744000114



ISO9001 Verification No.:
00217Q26901R1S



Guangzhou Ostec Electronic Technology Co., Limited
Manufacturer: No.8 West Lane, Jiangcheng Road, Bangjiang EastVillage, Shiqi Town,
Panyu District, Guangzhou, China.

www.ostec.com.cn

